

## VENAFLON HF TUBE

HOSES › PIPES suitable for ATEX environment

Chemical-resistant silicone hose with PFA inner layer compatible with highly aggressive chemicals.

### Applications

It is particularly suitable for the transport of liquids or semi-liquids chemically aggressive by impulse or suction in industry food, cosmetics, chemicals and pharmaceuticals.

It has a wide field of application thanks to its construction, which provides a balance between strength and lightness.


### Property

- Odorless, tasteless and completely non-toxic.
- High flexibility.
- White and smooth appearance of the inner PFA layer, translucent and smooth appearance of the outer silicone layer.
- Can be fitted with 316L stainless steel fittings on each ends with a roughness value lower than 0.8  $\mu\text{m}$  (or 0.5  $\mu\text{m}$  on request).
- On request it can be equipped with Clamp fittings with internal PFA covering.
- Operating temperature range  $-30^{\circ}\text{C}$  ( $-22^{\circ}\text{F}$ ) to  $+150^{\circ}\text{C}$  ( $302^{\circ}\text{F}$ ).
- The hose is produced with a maximum length of 20 m (65.62 ft).
- The vacuum resistance is 0.9 bar (13.05 psi).

**Construction This reference is manufactured with a white inner layer of PFA**

**(perfluoroalkoxy), the reinforcements are made of polyester and the spiral is made of stainless steel**

**(AISI 304) coated with silicone.**

**environment version  with conductive PFA inner layer, black color certified for food and pharmaceutical use (see below).**


### Characteristics

- Outer layer in VMQ Silicone, smooth and translucent.
- Reinforcement Stainless steel wire spring and polyester fabric reinforcement.
- Inner layer PFA fluoropolymer, smooth and white.
- Temperature  $-30^{\circ}\text{C}$   $+150^{\circ}\text{C}$
- Production length 20m. Can be cut to size upon request.
- Vacuum pressure 0.9 bar (13.05 psi).

### Regulations

The inner layer of PFA complies with :

- U.S. FDA Standard 21 CFR 177.1550
- USP Class VI in vivo test - USP Class VI in vitro test
- ISO 10993-5, 10 and 11 - Reg 1935/2004 / EEC and Reg 10/2011 / EEC



[Mastertubi.it/q?1035](http://Mastertubi.it/q?1035)

The outer silicone layer complies with :

- US FDA Standard 21 CFR 177.2600
- USP Class VI in vivo testing
- USP Class VI in vitro test
- ISO 10993-5, 10 and 11
- European Pharmacopoeia 3.1.9

<b>DIAMETER INTERNAL mm</b>	<b>WALL THICKNESS ISO 1307 +/- 0.80 mm</b>	<b>OPERATING PRESSURE ISO 1402 Bar at 20°C</b>	<b>RADIUS OF CURVATURE mm</b>
10	6	10	40
13	6	10	45
16	6	10	55
19	6	10	65
25	6	10	85
32	6	10	120
38	6.5	10	140
51	8	10	180
63.5	8	5	320
76	8	5	380
100	9	3	500